

Directions for Use



	<p>Align the Handle on the Hi-Flow Canister in between the alignment Pins.</p>
	<p>Remove the Tandem Port Cap as well as the large Pour Spout Cap on the Hi-Flow Canister cover.</p> <p>IMPORTANT: Be sure the caps are firmly affixed to the Patient Port and Vacuum Port.</p>
	<p>Rotate the Drainhead in the counterclockwise direction onto the top of the Hi-Flow Canister to seal drainhead to the Tandem and Pour ports. The drainhead will lock once it has been lowered into the correct position.</p>
	<p>Grasp the Handle and rotate the swingarm in the clockwise direction until the swingarm comes in contact with the stop. The canister will automatically drain.</p>
	<p>If you would like to rinse the interior of the canister, press the Rinse Button at this time.</p>
	<p>After the canister has drained, grasp the Handle and rotate the swingarm in the counterclockwise direction until the swingarm is in the home position.</p>
	<p>Grasp the Release Lever on the bottom of the drainhead and rotate the Drainhead clockwise until it rests against its stop. Recap all ports and dispose of according to hospital policy.</p>

Daily Maintenance - Run 3000cc of water/PH neutral enzymatic cleaner mixture through the unit at the end of the day.

For Service Contact:

Name: _____

Phone: _____