
SPECIMEN ADAPTOR

SPECIMEN ADAPTOR	4	ADATTATORE PER CAMPIONI	12
АДАПТЕР ЗА ПРОБИ	4	SPECIMEN ADAPTOR	12
PRILAGODNIK ZA UZORKE	5	PARAUGU ADAPTERS	13
ADAPTÉR NA VZORKY	6	MÈGINIO ADAPTERIS	14
PRØVEADAPTER	6	ADAPTER TAL-KAMPJUN	14
SPECIMENADAPTER	7	ADAPTER DO PRÓBEK	15
NÄIDISADAPTER	8	ADAPTADOR DE AMOSTRAS	16
NÄYTEADAPTERI	8	ADAPTOR PROBĀ	16
ADAPTATEUR D'ÉCHANTILLON	9	ADAPTÉR VZORKY	17
PROBENADAPTER	10	ADAPTER ZA VZORCE	18
ΠΡΟΣΑΡΜΟΓΕΑΣ ΔΕΙΓΜΑΤΟΣ	10	ADAPTADOR DE MUESTRAS	18
MINTAADAPTER	11	ADAPTER FÖR PROVBEHÅLLARE	19


SPECIMEN ADAPTOR

	EN	BG	HR	CZ	DA	NL
1	SOCK	РЪКАВ	ŽLIJEB ZA ISTAKANJE	TAMPON	SOK	SOK
2	PATIENT PORT	ПОРТ НА ПАЦИЕНТА	PRIKLJUČAK ZA PACIJENTA	PORT PACIENTA	PATIENTPORT	PATIËNTPOORT
3	CAP	ЗАПУШАЛКА	POKLOPAC	VÍČKO	HÆTTE	DOP


	ET	FI	FR	DE	EL	HU
1	PROOVIKOTT	SUOJUS	CHAUSSETTE	PROBENAUFFANGBEUTEL	ΚΑΛΥΜΜΑ	FOGLALAT
2	PATSIENDIAVA	POTILASKANAVA	PORT PATIENT	PATIENTENANSCHLUSS	ΘΥΡΑ ΑΣΘΕΝΗ	BETEGOLDALI NYÍLÁS
3	KORK	TULPPA	CAPUCHON	KAPPE	ΠΩΜΑ	KUPAK

	IT	GA	LV	LT	MT	PL
1	MANICHETTA	SOCK	IELIKTNIS	MAIŠELIS	KALZETTA	OKŁADZINA
2	PORTA PER PAZIENTE	PATIENT PORT	PACIENTA PIESLĒGVIETA	PACIENTO ANGA	PORT TAL-PAZJENT	PORT PACJENTA
3	CAPPUCCIO	CAP	VĀKS	DANGTELIS	KAPPA	ZĄSŁEPKA

	PT	RO	SK	SL	ES	SV
1	MANGA	MANȘON	VRECKO	PREVLEKA	CALCETÍN	SOCKA
2	PORTA DO DOENTE	PORT PENTRU PACIENT	PORT PACIENTA	ODPRTINA ZA PACIENTA	PUERTO DEL PACIENTE	PATIENTPORT
3	TAMPA	CAPAC	UZÁVER	POKROV	TAPA	LOCK

HI-FLOW COVER

HYDROPHOBIC COVER

	EN	BG	HR	CZ	DA	NL
1	PATIENT	ПАЦИЕНТ	PACIJENT	PACIENT	PATIENT	PATIËNT
2	VACUUM	ВАКУУМ	VAKUUM	PODTLAK	VAKUUM	VACUÛM
3	TANDEM	ДВОЕН	TANDEM	TANDEM	TANDEM	TANDEM
4	POUR SPOUT	ДОЗАТОР	NASTAVAK ZA IZLIJEVANJE	VYLÉVACÍ OTVOR	TAPSNUDE	SCHENKTUIT

	ET	FI	FR	DE	EL	HU
1	PATSIENT	POTILAS	PATIENT	PATIENT	ΑΣΘΕΝΗΣ	BETEG
2	VAAKUM	VAKUUMI	VIDE	VAKUUM	KENO	VÁKUUM
3	TANDEMAVA	KAKSOIS	TANDEM	TANDEM	ΠΑΡΑΛΛΗΛΗ ΘΥΡΑ	TANDEM
4	VALAMISTILA	KAATONOKKA	BECVERSEUR	AUSGIESSER	ΣΤΟΜΙΟ ΠΟΗΣ	KIÖNTÖNYÍLÁS

	IT	GA	LV	LT	MT	PL
1	PAZIENTE	PATIENT	PACIENTS	PACIENTAS	PAZJENT	PACJENT
2	VUOTO	VACUUM	VAKUUMS	VAKUUMAS	VAKWU	PODCIŚNIENIE
3	TANDEM	TANDEM	TANĒMS	TANDEMAS	TANDEM	TANDEM
4	BECCUCCIO DOSATORE	POUR SPOUT	LEJAMTEKNE	IŠPYLIMO ANGA	ŻENNUNA GHATTFERRIGH	WLEW

	PT	RO	SK	SL	ES	SV
1	DOENTE	PACIENT	PACIENT	PACIENT	PACIENTE	PATIENT
2	VÁCUO	VACUUM	VÁKUUM	VAKUUM	VACÍO	VAKUUM
3	TANDEM	TANDEM	TANDEM	TANDEM	TÁNDEM	TANDEM
4	BOCAL	PORT DE GOLIRE	DÁVKOVACIA HUBICA	DULEC ZA IZLIVANJE	BOQUILLA DE VERTIDO	HÄLLTAPP

SPECIMEN ADAPTOR

INTENDED USE: Intended for connection to suction canisters to allow for the collection of specimens removed during surgical or medical procedures.

DIRECTIONS FOR USE - ASSEMBLY

Note: Begin directions only after properly assembling suction canister.

1. Insert Sock into Pour Spout.
2. Inspect that Specimen Adaptor is sealed to Canister Cover.

3. Inspect all connections.
4. Turn on vacuum source and test for any leakage.
5. Collect desired sample.

REMOVING SAMPLE

1. Remove tubing from Specimen Adaptor and cap port.
2. Remove Specimen Adaptor using caution.

DISPOSAL

- Dispose of according to your facility's standard operating procedures for biohazard disposal and both federal and local disposal regulations.

TROUBLESHOOTING

- Loss of Vacuum: Check that vacuum is on, canister is properly sealed, and that all connections are tight and tubing is not kinked. If loss continues, replace canister.

⚠ CAUTION

- Intended User: To be used by professionally trained healthcare personnel.
- **⊗ Single Use Only. Do not attempt to clean, sterilize or reuse specimen adaptor.** Possible consequence of reuse include exposure to pathogens.
- **Do not** exceed vacuum level of 25" (640mm) Hg (84.7 kPa).
- **☠** Contents are considered potentially hazardous. Use appropriate PPE and handle accordingly.
- **☀** Store in a dark place. Long term exposure to light may compromise product performance and result in breakage during use.

Any serious incident related to the use of this product, should be reported to both the manufacturer and the health authority/competent authority where the product is established.

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REMOVING SAMPLE

1. Remove tubing from Specimen Adaptor and cap port.

2. Remove Specimen Adaptor using caution.

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⚠ CAUTION

- Intended User: To be used by professionally trained healthcare personnel.
-  **Single Use Only. Do not attempt to clean, sterilize or reuse specimen adaptor.** Possible consequence of reuse include exposure to pathogens.
- **Do not** exceed vacuum level of 25" (640mm) Hg (84.7 kPa).

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-  Contents are considered potentially hazardous. Use appropriate PPE and handle accordingly.
-  Store in a dark place. Long term exposure to light may compromise product performance and result in breakage during use.

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UK CA Responsible Person: Bemis Limited