



Control Today. Protect Tomorrow.™

**BEMIS**  
**Sharps Sentinel**  
Sharps Disposal System



**Sentinel**™  
by BEMIS®

WASTE CONTAINER CATALOG

## OUR COMMITMENT TO SUSTAINABLE PRODUCTS AND MANUFACTURING PROCESSES

Protecting the health of medical professionals, patients, medical waste handlers, communities and the environment is at the core of Bemis' mission.

Our commitment to protecting the environment extends not just to the products we bring to market but to our manufacturing processes as well. We work to control our carbon footprint through environmentally responsible operational practices, aiming to reduce eco-strain by minimizing impact from areas such as raw material consumption, production and disposal. We look to reprocessed materials whenever possible, and continue to reduce waste, energy, resource consumption and emissions during production. We also monitor our own progress through our Environmental Impact Committee to make sure we are moving in the right direction.

Through progress and innovation, Bemis Manufacturing Company will continue to be an environmentally-friendly neighbor and a responsible global citizen.

Learn more about environmental stewardship: [BemisHealthCare.com/customer-service](https://www.bemishealthcare.com/customer-service)



# FOR MORE THAN 50 YEARS, BEMIS HEALTH CARE HAS PROVIDED HIGH QUALITY, DISPOSABLE MEDICAL PRODUCTS FOR THE HEALTH CARE INDUSTRY

The Sentinel by Bemis™ waste container portfolio is available in a variety of sizes and clinical use types to meet biohazardous sharps and specialty waste container disposal that meet or exceed OSHA, NIOSH, CDC and JCAHO requirements. Units are puncture resistant, meeting ISO:23907-1 needle puncture standard.



## SHARPS CONTAINERS (BIOHAZARDOUS WASTE)

Sharps containers are designed with quality, puncture resistant materials that offer the safeguards health care facilities and clinicians deserve. Biohazard logo clearly distinguishes the container use from other disposal collection types. Available in multiple sizes, colors and lid openings to meet clinical use needs.

---



## CHEMOTHERAPY CONTAINERS

Chemotherapy containers are designed for safe collection and containment of trace chemotherapy waste. Trace waste is defined by CDC/NIOSH as "containing less than 3% by weight of original quantity of hazardous drug". Each container includes gasketed lids with absorbent pads. Available in three sizes with a choice of yellow or white containers.

---



## PHARMACEUTICAL CONTAINERS

Designed for the disposal of expired or unused non-hazardous pharmaceuticals. Each container includes gasketed lids with absorbent pads. Available in three blue container sizes with labeling to easily distinguish between other disposal collection types.

---



## HAZARDOUS PHARMACEUTICAL CONTAINERS

Hazardous pharmaceutical containers, commonly known as RCRA (Resource Conservation and Recovery Act) are designed to meet EPA requirements of P and U-listed chemical disposal. Each container includes gasketed lids with absorbent pads. The containers are available in three black container sizes with labeling to easily distinguish between other disposal collection types.

# BEMIS HEALTH CARE WASTE DISPOSAL SOLUTIONS

DISPOSAL TYPE	USE TYPE	CAPACITY	SAMPLE WASTE	LID TYPE	CATALOG #	COLOR
<b>SHARPS* (BIOHAZARD)</b>	Phlebotomy Sharps	1-Quart	1-10 ml syringes, blood collection tubes, pipettes	Dual openings/vertical drop	100 030	 Red
	Patient/Exam Room Sharps	5-Quart	1-60 ml syringes with needle	Rotating cylinder/horizontal drop	175 030	 Red
					175 020	 Beige
	Multi-Purpose Sharps	2-Gallon	1-60 ml syringes with needle, blood collection tubes, slides, pipettes	Angled chute/horizontal drop	102 030	 Red
					102 020	 Beige
	Patient/Exam Room Sharps	3-Gallon	1-60 ml syringes with needle	Rotating cylinder/horizontal drop	333 030	 Red
					333 020	 Beige
Multi-Purpose Sharps	3 Gallon	1-60 ml syringes with needle, blood collection tubes, slides, pipettes	Angled chute/horizontal drop	303 030	 Red	
				303 020	 Beige	
Multi-Purpose/ Large Volume Sharps	8-Gallon	Various sharps and larger OPIM items	Hinged for large items/separate opening for smaller items	108 030	 Red	
	11-Gallon			111 030	 Red	
<b>TRACE CHEMOTHERAPY</b>	Chemotherapy Medical Waste	2-Gallon	Various trace chemo sharps, vials, syringes	Round screw-top with 4 inch opening	202 004	 White
					202 040	 Yellow
	Chemotherapy Large Volume Medical Waste	8-Gallon	Various trace chemo sharps, tubing, vials, syringes, pads, gowns	Gasketed hinge for large items/separate opening for smaller items	208 004	 White
					208 040	 Yellow
					211 004	 White
	11-Gallon			211 040	 Yellow	
<b>NON-HAZARDOUS PHARMACEUTICAL</b>	Pharmaceutical Medical Waste	2-Gallon	NO SHARPS Non-hazardous pharmaceutical syringes/ vials with residual medication, small items	Round screw-top with 4 inch opening	4002 050	 Blue
	Pharmaceutical Large Volume Medical Waste	8-Gallon	NO SHARPS Non-hazardous pharmaceutical syringes/ vials with residual medication, tablets, creams, eyedrops, etc.	Hinged for large items/separate opening for smaller items	4008 050	 Blue
		11-Gallon			4011 050	 Blue
<b>HAZARDOUS PHARMACEUTICAL</b>	Hazardous Pharmaceutical Medical Waste	2-Gallon	NO SHARPS Small hazardous pharmaceutical P and U-listed items	Round screw-top with 4 inch opening	5002 070	 Black
	Hazardous Pharmaceutical Medical Waste	8-Gallon	NO SHARPS Small hazardous pharmaceutical P and U-listed items - inhalers, nitroglycerine	Hinges for large items/separate opening for smaller items	5008 070	 Black
		11-Gallon			5011 070	 Black

Yellow biohazard sharps containers for Canadian market not listed in table - refer to pages 6-8 for additional information

# GENERAL MEDICAL WASTE DISPOSAL & SEGREGATION GUIDE

Check your state guidelines for specific medical waste disposal regulations

		
<p><b>NON-HAZARDOUS LANDFILL</b></p>	<p><b>AUTOCLAVE/INCINERATE</b></p>	
<p><b>NON-HAZARDOUS BAG WASTE</b></p>	<p><b>BIOHAZARDOUS BAG WASTE</b></p>	<p><b>BIOHAZARDOUS (SHARPS)</b></p>
<p>IV bags and tubing empty med Vials trash/wrappers Dressings chux Empty drainage bags Disposable items Sanitary napkins Batteries</p>	<p>Infectious materials (OPIM) Blood tubing, bags, hemovacs, pleurevacs Blood soaked or dripping dressings Contaminated waste from isolation patients Suction canisters and liners with bodily fluids</p>	<p>Needles Glass slides, pipettes Scalpels Bloody broken glass Used sutures Empty syringes Trocars, introducers, catheters, electrocautery tips</p>
<p><b>WHITE OR CLEAR BAGS</b></p>	<p><b>RED BAG WITH BIOHAZARD SYMBOL</b></p>	<p><b>CONTAINER WITH BIOHAZARD SYMBOL</b></p>
		
<p><b>LOW TEMP INCINERATION</b></p>	<p><b>HIGH TEMP INCINERATION</b></p>	
<p><b>CHEMOTHERAPY (TRACE)</b></p>	<p><b>NON-HAZARDOUS PHARMACEUTICAL</b></p>	<p><b>HAZARDOUS PHARMACEUTICAL</b></p>
<p>Trace chemo waste of no more than 3% weight of original components All supplies used to make and administer chemo medication including sharps Examples: Tubing, empty bags/ bottles/vials, syringes, needles, gloves, pads masks, gowns, wipes, etc.</p>	<p>NO SHARPS Residual (pourable) meds Discontinued IV bags/tubing Partially used prescription or OTC medication Examples: Vials, tablets, capsules, powders, contained liquids, creams/ lotions, eyedrops, residual or wasted narcotics, unused narcotic patches</p>	<p>NO SHARPS RCRA Hazardous P and U-listed chemicals Examples: Inhalers, nitroglycerin tablets, unused acetone, bulk chemotherapy waste</p>
<p><b>CYTOTOXIC &amp; BIOHAZARD SYMBOLS</b></p>	<p><b>NON-HAZARDOUS PHARMACEUTICAL LABEL</b></p>	<p><b>HAZARDOUS PHARMACEUTICAL LABEL</b></p>

# PATIENT & EXAM ROOM SHARPS CONTAINERS

Designed for use in patient rooms, physicians' offices, clinics and other high-traffic areas



## FEATURES & BENEFITS

- Rotating lid design restricts access to container contents
- Lid design accommodates syringes up to 60ccs
- Horizontal drop allows maximum capacity of sharps
- Translucent containers and lids help identify fill level
- Can be used on tabletop or with locking cabinets or bracket wall mounting options for safety and convenience
- Containers lock for final disposal
- Nestable to save space and reduce shipping costs
- English, Spanish and French labeling

*Note: For final disposal - render non-infectious by autoclave then incinerate or grind/shred*

## PATIENT & EXAM ROOM SHARPS CONTAINERS WITH MOUNTING OPTIONS

### PRODUCT/DESCRIPTION

### CATALOG #

### PACKAGING

### MOUNTING OPTIONS



#### 5-QUART SHARPS CONTAINER

- rotating cylinder top that allows for hands-free disposal
- built-in handles for easy transport
- uses bracket/key #445 or cabinet 1051 for wall mounting

Height: 10" Length: 11" Width: 5 1/4"

- 175 020 Translucent Beige
- 175 030 Translucent Red
- 175 040\* Translucent Yellow

32/CASE



1051



445 020



#### 3-GALLON SHARPS CONTAINER

- rotating cylinder top that allows for hands-free disposal
- built-in handles for easy transport
- uses bracket #423 or 425, or cabinet 1031 for wall mounting

Height: 15" Length: 13 7/8" Width: 6 7/8"

- 333 020 Translucent Beige
- 333 030 Translucent Red
- 333 040\* Translucent Yellow

12/CASE



1031



423 000



425 000

# PHLEBOTOMY & MULTI-PURPOSE SHARPS CONTAINERS

Designed for use in restricted areas such as nurses' stations, med carts, laboratories, operating rooms and pharmacies



## FEATURES & BENEFITS

Available in multiple volume sizes, from 1-quart to 11-gallon to accommodate clinical needs

Built-in handles for easy transport (excluding phlebotomy)

Translucent containers and lids help identify fill level on 1-quart, 2 and 3-gallon sizes

2 and 3-gallon containers can be wall-mounted or used freestanding on tabletop or cart

8 and 11-gallon have large lid openings to accommodate a variety of biohazard sharps and materials and available with rolling cart for mobility and transport

Lids can be closed for temporary or final transport closure

All containers nestable for space savings

Horizontal drop on most units maximizes container volume

English, Spanish and French labeling

*Note: For final disposal - render non-infectious by autoclave then incinerate or grind/shred*

## PHLEBOTOMY & MULTI-PURPOSE SHARPS CONTAINERS WITH MOUNTING OPTIONS

PRODUCT/DESCRIPTION	CATALOG #	PACKAGING	MOUNTING OPTIONS
 <p><b>1-QUART PHLEBOTOMY CONTAINER</b></p> <ul style="list-style-type: none"> <li>fits phlebotomy trays</li> <li>dual openings, vertical drop</li> <li>uses holder #400 to secure container on tray or tabletop</li> </ul> <p>Height: 6 7/16" Length: 4 1/4" Width: 4 1/4"</p>	<ul style="list-style-type: none"> <li>100 030 Translucent Red</li> <li>100 040* Translucent Yellow</li> </ul>	100/CASE	 <p>400 000</p>
 <p><b>2-GALLON SHARPS CONTAINER</b></p> <ul style="list-style-type: none"> <li>accommodates syringes up to 60cc</li> <li>built-in handles for easy transport</li> <li>uses bracket/key #445 for wall mounting</li> </ul> <p>Height: 8 5/8" Length: 11 5/8" Width: 7 3/4"</p>	<ul style="list-style-type: none"> <li>102 020 Translucent Beige</li> <li>102 030 Translucent Red</li> <li>102 040* Translucent Yellow</li> </ul>	30/CASE	 <p>445 020</p>

\* Available for Canadian market

PHLEBOTOMY & MULTI-PURPOSE SHARPS CONTAINERS WITH MOUNTING OPTIONS CONTINUED

PRODUCT/DESCRIPTION

CATALOG #

PACKAGING

MOUNTING  
OPTIONS



**3-GALLON SHARPS CONTAINER**

- accommodates syringes up to 60cc
- built-in handles for easy transport
- uses bracket #423 or 425 for wall mounting

**Height:** 13 ½" **Length:** 13 ⅞" **Width:** 6 ⅞"

■ 303 020  
Translucent Beige


■ 303 030  
Translucent Red

■ 303 040\*  
Translucent Yellow

12/CASE



423 000



425 000



**8-GALLON SHARPS CONTAINER**

- hinged lid accommodates large items
- separate round opening accepts small items
- molded-in handles for easy lifting and transport
- uses cart #430 for mobility

**Height:** 15 ⅞" **Length:** 16 ½" **Width:** 11 ¼"

■ 108 030  
Red

■ 108 040\*  
Yellow

10/CASE



430 030



**11-GALLON SHARPS CONTAINER**

- hinged lid accommodates large items
- separate round opening accepts small items
- molded-in handles for easy lifting and transport
- uses cart #430 for mobility

**Height:** 22 ½" **Length:** 16 ½" **Width:** 11 ¼"

■ 111 030  
Red

■ 111 040\*  
Yellow

6/CASE



430 030

\* Available for Canadian market

# CHEMOTHERAPY SHARPS CONTAINERS

Designed for safe collection and containment of trace chemotherapy waste



## FEATURES & BENEFITS

Designed for safe collection and containment of trace chemotherapy waste including needles, syringes, IVs, empty or partially empty vials, as well as gloves and gowns

Gasketed lids and absorbent spill pads provide leak resistance

Available in white or yellow colors in 2, 8 and 11-gallon sizes

Clearly labeled with cytotoxic chemotherapy warnings and biohazard symbols

Nestable to save space and reduce shipping costs

English, Spanish and French labeling

*Note: Chemotherapy containers must be incinerated*

## CHEMOTHERAPY SHARPS CONTAINERS WITH MOUNTING OPTIONS

PRODUCT/DESCRIPTION	CATALOG #	PACKAGING	MOUNTING OPTIONS
 <p><b>2-GALLON CHEMOTHERAPY CONTAINER</b></p> <ul style="list-style-type: none"> <li>• screw-on cap; posts on lid for holding cap in place when in use</li> <li>• for use on tabletop or wall-mounted</li> <li>• uses bracket/key #445 for wall mounting</li> </ul> <p><b>Height: 9" Length: 11 5/8" Width: 7 3/4"</b></p>	<p><input type="checkbox"/> 202 004 White w/gasket lid/ absorbent pad</p> <p><input checked="" type="checkbox"/> 202 040 Yellow w/gasket lid/ absorbent pad</p>	30/CASE	 <p>445 020</p>
 <p><b>8-GALLON CHEMOTHERAPY CONTAINER</b></p> <ul style="list-style-type: none"> <li>• hinged lid accommodates large items</li> <li>• separate round opening accepts small items</li> <li>• molded-in handles for easy lifting and transport</li> <li>• uses cart #430 for mobility</li> </ul> <p><b>Height: 15 7/8" Length: 16 1/2" Width: 11 13/16"</b></p>	<p><input type="checkbox"/> 208 004 White w/gasket lid/ absorbent pad</p> <p><input checked="" type="checkbox"/> 208 040 Yellow w/gasket lid/ absorbent pad</p>	10/CASE	 <p>430 030</p>
 <p><b>11-GALLON CHEMOTHERAPY CONTAINER</b></p> <ul style="list-style-type: none"> <li>• hinged lid accommodates large items</li> <li>• separate round opening accepts small items</li> <li>• molded-in handles for easy lifting and transport</li> <li>• uses cart #430 for mobility</li> </ul> <p><b>Height: 22 1/2" Length: 16 1/2" Width: 11 13/16"</b></p>	<p><input type="checkbox"/> 211 004 White w/gasket lid/ absorbent pad</p> <p><input checked="" type="checkbox"/> 211 040 Yellow w/gasket lid/ absorbent pad</p>	6/CASE	 <p>430 030</p>

# PHARMACEUTICAL WASTE CONTAINERS

Designed for the disposal of expired or unused non-hazardous pharmaceuticals



## FEATURES & BENEFITS

Designed for the disposal of expired or unused nonhazardous pharmaceuticals such as antibiotics, IVs, pills, etc.

Ideal for use in pharmacies, nursing stations and med carts

Blue color and special labeling easily distinguishes from other disposal containers for sharps, chemotherapy and Hazardous pharmaceutical waste

Gasketed lids and absorbent pads provide leak resistance

Available in 2, 8 and 11-gallon sizes

Nestable to save space and reduce shipping costs

English, Spanish and French labeling

*Notes: Container **NOT** for sharps use; pharmaceutical waste must be incinerated*

## PHARMACEUTICAL WASTE CONTAINERS WITH MOUNTING OPTIONS

PRODUCT/DESCRIPTION	CATALOG #	PACKAGING	MOUNTING OPTIONS
 <p><b>2-GALLON PHARMACEUTICAL CONTAINER</b></p> <ul style="list-style-type: none"> <li>• screw-on cap; posts on lid for holding cap in place when in use</li> <li>• for use on tabletop or wall-mounted</li> <li>• uses bracket/key #445 for wall mounting</li> </ul> <p><b>Height: 9" Length: 11 5/8" Width: 7 3/4"</b></p>	<p>■ 4002 050 blue w/gasket lid/ absorbent pad</p>	<p>30/CASE</p>	 <p>445 020</p>
 <p><b>8-GALLON PHARMACEUTICAL CONTAINER</b></p> <ul style="list-style-type: none"> <li>• hinged lid accommodates large items</li> <li>• separate round opening accepts small items</li> <li>• molded-in handles for easy lifting and transport</li> <li>• uses cart #430 for mobility</li> </ul> <p><b>Height: 15 7/8" Length: 16 1/2" Width: 11 13/16"</b></p>	<p>■ 4008 050 blue w/gasket lid/ absorbent pad</p>	<p>10/CASE</p>	 <p>430 030</p>
 <p><b>11-GALLON PHARMACEUTICAL CONTAINER</b></p> <ul style="list-style-type: none"> <li>• hinged lid accommodates large items</li> <li>• separate round opening accepts small items</li> <li>• molded-in handles for easy lifting and transport</li> <li>• uses cart #430 for mobility</li> </ul> <p><b>Height: 22 1/2" Length: 16 1/2" Width: 11 13/16"</b></p>	<p>■ 4011 050 blue w/gasket lid/ absorbent pad</p>	<p>6/CASE</p>	 <p>430 030</p>

# HAZARDOUS PHARMACEUTICAL WASTE CONTAINERS

Designed for the disposal of hazardous pharmaceuticals



## FEATURES & BENEFITS

- Container designed to meet EPA requirements of P and U-listed chemicals
- Gasketed lids and absorbent pads provide leak resistance
- Lids can be closed for temporary or final transport closure
- Available in 2, 8 and 11-gallon sizes
- Nestable to save space and reduce shipping costs
- English, Spanish and French labeling

*Notes: Container **NOT** for sharps use; RCRA waste must be incinerated by a federally permitted hazardous waste incinerator*

## HAZARDOUS PHARMACEUTICAL WASTE CONTAINERS WITH MOUNTING OPTIONS

PRODUCT/DESCRIPTION	CATALOG #	PACKAGING	MOUNTING OPTIONS
 <p><b>2-GALLON RCRA CONTAINER</b></p> <ul style="list-style-type: none"> <li>screw-on cap; posts on lid for holding cap in place when in use</li> <li>for use on tabletop or wall-mounted</li> <li>uses bracket/key #445 for wall mounting</li> </ul> <p><b>Height: 9" Length: 11 5/8" Width: 7 3/4"</b></p>	<p>■ 5002 070 black w/gasket lid/ absorbent pad</p>	<p>30/CASE</p>	 <p>445 020</p>
 <p><b>8-GALLON RCRA CONTAINER</b></p> <ul style="list-style-type: none"> <li>hinged lid accommodates large items</li> <li>separate round opening accepts small items</li> <li>molded-in handles for easy lifting and transport</li> <li>uses cart #430 for mobility</li> </ul> <p><b>Height: 15 7/8" Length: 16 1/2" Width: 11 13/16"</b></p>	<p>■ 5008 070 black w/gasket lid/ absorbent pad</p>	<p>10/CASE</p>	 <p>430 030</p>
 <p><b>11-GALLON RCRA CONTAINER</b></p> <ul style="list-style-type: none"> <li>hinged lid accommodates large items</li> <li>separate round opening accepts small items</li> <li>molded-in handles for easy lifting and transport</li> <li>uses cart #430 for mobility</li> </ul> <p><b>Height: 22 1/2" Length: 16 1/2" Width: 11 13/16"</b></p>	<p>■ 5011 070 black w/gasket lid/ absorbent pad</p>	<p>6/CASE</p>	 <p>430 030</p>

# SHARPS CABINETS, WALL MOUNTS & ACCESSORIES



## CABINETS & ACCESSORIES

PRODUCT/DESCRIPTION	CATALOG #	PACKAGING
 <p><b>5-QUART WALL CABINET</b></p> <ul style="list-style-type: none"> <li>• two universal keys provided</li> <li>• cabinet back accommodates most competitive drilling/wall-mounting pattern</li> <li>• for use with sharps container #175</li> </ul> <p><b>Height: 11 ½" Length: 13" Width: 6"</b></p>	 <p>1051</p>	1/CASE
 <p><b>3-GALLON WALL CABINET</b></p> <ul style="list-style-type: none"> <li>• two universal keys provided</li> <li>• cabinet back accommodates most competitive drilling/wall-mounting pattern</li> <li>• for use with sharps container #333</li> </ul> <p><b>Height: 16 ¾" Length: 15" Width: 7"</b></p>	 <p>1031</p>	1/CASE
 <p><b>GLOVE BOX HOLDER</b></p> <ul style="list-style-type: none"> <li>• snaps easily onto bottom of cabinets</li> <li>• for use with wall cabinets #1051 and 1031</li> </ul> <p><b>Height: 5 ½" Length: 11 ⅝" Width: 4 ⅜"</b></p>	 <p>1041</p>	1/CASE

## BRACKETS & ACCESSORIES

 <p><b>MOUNTING BRACKET/KEY FOR 5-QUART &amp; 2-GALLON</b></p> <ul style="list-style-type: none"> <li>• mounting hardware included with 10 keys</li> <li>• compatible with all 5 Quart and 2 Gallon containers to hold securely to wall</li> </ul> <p><b>Height: 6" Length: 10 ¼" Width: 2 ¼"</b></p>	445 020	10/CASE
 <p><b>MOUNTING BRACKET FOR 5-QUART &amp; 2-GALLON</b></p> <ul style="list-style-type: none"> <li>• mounting hardware included</li> <li>• compatible with all 5 Quart and 2 Gallon containers to hold securely to wall</li> <li>• requires key #440 to remove container from bracket</li> </ul> <p><b>Height: 6" Length: 10 ¼" Width: 2 ¼"</b></p>	435 020	10/CASE
 <p><b>KEY FOR #435 MOUNTING BRACKET</b></p> <ul style="list-style-type: none"> <li>• plastic key removes #175 5-quart and all 2-gallon containers from bracket</li> <li>• for use with mounting bracket #435</li> </ul> <p><b>Height: 3 ⅜" Length: 1" Width: ⅝"</b></p>	440 020	10/CASE

## BRACKETS & ACCESSORIES CONTINUED

PRODUCT/DESCRIPTION	CATALOG #	PACKAGING
 <p><b>NON-LOCKING WIRE BRACKET FOR 3-GALLON</b></p> <ul style="list-style-type: none"> <li>holds sharps containers #303 and 333 securely to wall</li> <li>mounting hardware included</li> </ul> <p><b>Height: 7 ½" Length: 12" Width: 6"</b></p>	423 000	5/CASE
 <p><b>LOCKING WIRE BRACKET FOR 3-GALLON</b></p> <ul style="list-style-type: none"> <li>holds sharps containers #303 and 333 securely to wall</li> <li>same as bracket #423 but with swing arm that fits over container and secures with padlock</li> <li>mounting hardware included</li> <li>requires padlock/key #436</li> </ul> <p><b>Height: 9 ½" Length: 12" Width: 6"</b></p>	425 000	5/CASE
 <p><b>PADLOCK/KEY</b></p> <ul style="list-style-type: none"> <li>for use with locking wire bracket #425</li> </ul> <p><b>Height: 1 ¾" Length: 1 ½" Width: ½"</b></p>	426 000	1/CASE

## CARTS & HOLDERS



**ROLLING CART FOR 8 & 11-GALLON CONTAINERS**

- durable construction
- accommodates all 8 and 11-gallon containers

**Height: 16" Length: 17" Width: 13"**

430 030

1/CASE



**PHLEBOTOMY HOLDER**

- translucent white allows visibility of container label
- adhesive on bottom secures container to tabletop or tray
- for use with phlebotomy container #100

**Height: 6 ⅝" Length: 4 ¼" Width: 4 ¼"**

400 000

20/CASE

# CONTAINER MOUNTING OPTIONS QUICK REFERENCE TABLE

CONTAINER CATALOG #	COLOR	SIZE	MOUNTING TYPE & CATALOG #						
			5-QT SHARPS CABINET	3-GAL SHARPS CABINET	PHLEBOTOMY HOLDER	WIRE BRACKET (NON-LOCKING)	WIRE BRACKET (LOCKING)*	WALL BRACKET & KEY**	ROLLING CART
			1051	1031	400 000	423 000	425 000	445 020	430 030
<b>SHARPS CONTAINERS</b>									
100 040	Yellow	1-Quart			•				
100 030	Red	1-Quart			•				
175 040	Yellow	5-Quart	•					•	
175 030	Red	5-Quart	•					•	
175 020	Beige	5-Quart	•					•	
102 040	Yellow	2-Gallon						•	
102 030	Red	2-Gallon						•	
102 020	Beige	2-Gallon						•	
303 040	Yellow	3-Gallon				•	•		
303 030	Red	3-Gallon				•	•		
303 020	Beige	3-Gallon				•	•		
333 040	Yellow	3-Gallon		•		•	•		
333 030	Red	3-Gallon		•		•	•		
333 020	Beige	3-Gallon		•		•	•		
108 030	Red	8-Gallon							•
108 040	Yellow	8-Gallon							•
111 030	Red	11-Gallon							•
111 040	Yellow	11-Gallon							•
<b>CHEMOTHERAPY CONTAINERS</b>									
202 004	White	2-Gallon						•	
202 040	Yellow	2-Gallon						•	
208 004	White	8-Gallon							•
208 040	Yellow	8-Gallon							•
211 004	White	11-Gallon							•
211 040	Yellow	11-Gallon							•
<b>NON-HAZARDOUS PHARMACEUTICAL CONTAINERS</b>									
4002 050	Blue	2-Gallon						•	
4008 050	Blue	8-Gallon							•
4011 050	Blue	11-Gallon							•
<b>HAZARDOUS PHARMACEUTICAL CONTAINERS (RCRA)</b>									
5002 070	Black	2-Gallon						•	
5008 070	Black	8-Gallon							•
5011 070	Black	11-Gallon							•

\* Purchase padlock and key separately - catalog # 426 000

\*\* Wall brackets and keys can be purchased separately - catalog # 435 020 and 440 020

# NOTES

A series of horizontal dotted lines for writing notes.

Products included in this catalog are available for use within the United States and Canadian market.

For more information on medical waste disposal & Bemis Health Care  
[BemisHealthCare.com](http://BemisHealthCare.com)

Waste regulations vary by state. For a complete list of state regulatory agencies and resources go to the following link within the EPA's website: [www.epa.gov/rcra/medical-waste](http://www.epa.gov/rcra/medical-waste)

Control Today. Protect Tomorrow.™

**Bemis Manufacturing Company**

Health Care Division | 300 Mill Street | P.O. Box 901 | Sheboygan Falls, Wisconsin 53085 USA

Email: [Customer.Service@BemisHealthCare.com](mailto:Customer.Service@BemisHealthCare.com) | Phone: 920.467.4621 | Toll-Free Phone: 800.558.7651 | Fax: 800.499.8160

3/2022 Rev 02

©2022 0B7752112

