



# QUICK-DRAIN™ | Disposal Systems for Liquid Infectious Waste

*Now with Large Volume Feature*

Reduce Costs



Safeguard Your Staff



Reduce Waste

- **LARGE VOLUME ADAPTER.** Provides an option for draining large volume canisters, kick buckets, or other large volume containers.
- **EASY OPERATION.** Just insert dip tube into large volume canister and turn unit on.
- **QUICK.** Drains standard Bemis canisters at a rate of 500cc per second; Drains any large volume container up to 200cc per second.
- **SAFE.** Satisfies OSHA's Bloodborne Pathogen Standard to "provide engineering controls that isolate or remove the BBP from the workplace".
- **COST EFFECTIVE.** Reduces red bag waste, eliminates solidifiers. Low maintenance, no electricity needed, connects to existing plumbing.

**Model QDHF3000**

Quick-Drain for Bemis Hi-Flow Hard Canisters, 1/Case.

**Model QDQF3000**

Quick-Drain for Bemis Quick-Fit Suction Liners, 1/Case.

**Model QDHFLV3500**

Quick-Drain Hi-Flow with Large Volume Adapter, 1/Case.

**Model QDQFLV3500**

Quick-Drain Quick-Fit with Large Volume Adapter, 1/Case.

**Model QDLVRH**

Reusable Drain Hose for Large Volume Units, 2/Case.

**Model QDLVDT15**

15" Disposable Dip Tube for Large Volume Units, 100/Case.

**Model QDLVDT19**

19" Disposable Dip Tube for Large Volume Units, 100/Case.