



BRACKETS & ACCESSORIES




400 000
PHLEBOTOMY HOLDER
20/CASE


- Translucent White allows visibility of container label
- Adhesive on bottom secures container to tabletop or tray
- Dimensions: 6 5/16" H x 4 1/4" L x 4 1/4" W
- *For use with phlebotomy container #100*


410 020
MOUNTING BRACKET FOR WALLSAFE® SHARPS CONTAINERS
10/CASE


- Holds WallSafes #125 and 150 securely to wall
- Concealed by container when in use
- Mounting hardware included
- Dimensions: 6" H x 9" L x 1/2" W
- *Requires key #415 to remove container from bracket*


415 020
KEY FOR #410 MOUNTING BRACKET
10/CASE

- Unique plastic key removes WallSafes from bracket
- Dimensions: 3 1/4" H x 3/4" L x 3/4" W
- *For use with mounting bracket #410*


423 000
NON-LOCKING WIRE BRACKET FOR 3-GALLON
5/CASE

- Holds sharps containers #303 and 333 securely to wall
- Mounting hardware included
- Dimensions: 7 1/2" H x 12" L x 6" W
- *Can be mounted separately or slides onto mounting bracket #410*


425 000
LOCKING WIRE BRACKET FOR 3-GALLON
5/CASE

- Holds sharps containers #303 and 333 securely to wall
- Same as bracket #423 but with swing arm that fits over container and secures with padlock
- Mounting hardware included
- Dimensions: 9 1/2" H x 12" L x 6" W
- *Can be mounted separately or slides onto mounting bracket #410*
- *Requires padlock/key #426*


426 000
PADLOCK & KEY
1/CASE

- Dimensions: 1 3/4" H x 1 1/8" L x 1/2" W
- *For use with locking wire bracket #425*


430 030
ROLLING CART FOR 8- & 11-GALLON CONTAINERS
1/CASE

- Durable casters
- Sturdy construction
- Dimensions: 16" H x 17" L x 13" W
- *Accommodates all 8- and 11-gallon containers*


435 020
MOUNTING BRACKET FOR 5-QUART & 2-GALLON
10/CASE

- Dimensions: 6" H x 10 1/4" L x 2 1/4" W
- Mounting hardware included
- *Holds sharps container #175 5-quart and all 2-gallon containers securely to wall*
- *Requires key #440 to remove container from bracket*


440 020
KEY FOR #435 MOUNTING BRACKET
10/CASE

- Unique plastic key removes #175 5-quart and all 2-gallon containers from bracket
- Dimensions: 3 1/8" H x 1" L x 5/8" W
- *For use with mounting bracket #435*