



FEATURES

- Compact design uses minimum floor space.
- Large opening lids accommodate a variety of sharps.
- Lids provide temporary and final closure.
- Nestable to save space and reduce shipping costs.
- Molded-in handles for easy lifting and transport.
- Multilingual labeling ensures user safety.

LARGE VOLUME SHARPS CONTAINERS

Designed for high-volume use in operating rooms, pharmacies, laboratories, emergency rooms, dialysis rooms, and other areas. A rolling cart is available for convenient mobility and transport.





■ 108 030
■ 108 040

8-GALLON SHARPS CONTAINER

10/CASE

Red

Yellow

- Separate round opening accepts small items
- Molded-in handles for easy lifting and transport
- Hinged lid accommodates large items
- Dimensions: 15 ⁷/₈" H x 16 ¹/₂" L x 11 ¹³/₁₆" W
- *Uses cart #430 for mobility*



■ 111 030
■ 111 040

11-GALLON SHARPS CONTAINER

6/CASE

Red

Yellow

- Separate round opening accepts small items
- Molded-in handles for easy lifting and transport
- Hinged lid accommodates large items
- Dimensions: 22 ¹/₂" H x 16 ¹/₂" L x 11 ¹³/₁₆" W
- *Uses cart #430 for mobility*