



## FEATURES

Rotating lid design prevents access to container contents.

Horizontal drop allows maximum capacity of sharps.

Translucent containers and lids help identify fill level.

Containers lock for final disposal.

Accommodate syringes up to 60cc.

Nestable to save space and reduce shipping costs.

Can be autoclaved or incinerated for final disposal.

Multilingual labeling ensures user safety.

## PATIENT ROOM & EXAM ROOM SHARPS CONTAINERS

Designed for use in patient rooms, physicians' offices, clinics and other high-traffic areas.





-  125 020
-  125 030
-  125 040

**3-QUART WALLSAFE® SHARPS CONTAINER** **24/CASE**

**Translucent Beige**  
**Translucent Red**  
**Translucent Yellow**

- Rotating cylinder top allows one-handed disposal
- Ideal for low-volume areas
- Dimensions: 8" H x 11" L x 4 1/4" W
- *Requires use of bracket/key #410 and 415 for wall mounting*



-  150 020
-  150 030
-  150 040

**5-QUART WALLSAFE® SHARPS CONTAINER** **24/CASE**

**Translucent Beige**  
**Translucent Red**  
**Translucent Yellow**

- Rotating cylinder top allows one-handed disposal
- Dimensions: 12" H x 11" L x 4 1/4" W
- *Requires use of bracket/key #410 and 415 for wall mounting*



-  175 020
-  175 030
-  175 040

**5-QUART SHARPS CONTAINER** **32/CASE**

**Translucent Beige**  
**Translucent Red**  
**Translucent Yellow**

- Rotating cylinder top allows hands-free disposal
- Built-in handles for easy transport
- Dimensions: 10" H x 11" L x 5 1/4" W
- *Uses bracket/key #435 and 440 or cabinet 1051 for wall mounting*



-  333 020
-  333 030
-  333 040

**3-GALLON SHARPS CONTAINER** **12/CASE**

**Translucent Beige**  
**Translucent Red**  
**Translucent Yellow**

- Rotating cylinder top allows hands-free disposal
- Built-in handles for easy transport
- Dimensions: 15" H x 13 7/8" L x 6 7/8" W
- *Uses bracket #423 or 425, or cabinet 1031 for wall mounting*