



FEATURES

Gasketed lids and absorbent pads provide leak resistance.

Available in three sizes.

Nestable to save space and reduce shipping costs.

NOT for sharps.

Multilingual labeling ensures user safety.

Note: Pharmaceutical waste must be incinerated.

PHARMACEUTICAL WASTE CONTAINERS

Designed for the disposal of expired or unused non-hazardous pharmaceuticals such as antibiotics, IVs, ibuprofen, etc. Ideal for use in pharmacies, nursing stations and med carts. Blue color and special labeling easily distinguishes from other disposal containers for sharps, chemotherapy and RCRA waste.





■ 4002 050

2-GALLON PHARMACEUTICAL CONTAINER
Blue w/gasket lid/absorbent pad

30/CASE

- Screw-on cap; Posts on lid for holding cap in place when in use
- For use on tabletop or wall-mounted
- Dimensions: 9" H x 11 ⁵/₈" L x 7 ³/₄" W
- *Uses bracket/key #435 and 440 for wall mounting*



■ 4008 050

8-GALLON PHARMACEUTICAL CONTAINER
Blue w/gasket lid/absorbent pads

10/CASE

- Separate round opening accepts small items
- Hinged lid accommodates large items
- Molded-in handles for easy lifting and transport
- Dimensions: 15 ⁷/₈" H x 16 ¹/₂" L x 11 ¹³/₁₆" W
- *Uses cart #430 for mobility*



■ 4011 050

11-GALLON PHARMACEUTICAL CONTAINER
Blue w/gasket lid/absorbent pads

6/CASE

- Separate round opening accepts small items
- Hinged lid accommodates large items
- Molded-in handles for easy lifting and transport
- Dimensions: 22 ¹/₂" H x 16 ¹/₂" L x 11 ¹³/₁₆" W
- *Uses cart #430 for mobility*