



## FEATURES

Accommodate a large variety of sharps sizes.

Temporary and final closure.

Built-in handles for easy transport (excluding phlebotomy).

Translucent containers and lids help identify fill level.

Nestable to save space and reduce shipping costs.

Horizontal drop on most units maximizes container volume.

Can be autoclaved or incinerated for final disposal.

Multilingual labeling ensures user safety.

## PHLEBOTOMY & MULTI-USE SHARPS CONTAINERS

Designed for use in restricted areas such as nurses' stations, med carts, laboratories, ORs and pharmacies. Can be wall-mounted or used freestanding on tabletop or cart.





 **100 030**  
 **100 040**

**1-QUART PHLEBOTOMY CONTAINER**

**100/CASE**

**Translucent Red**

**Translucent Yellow**

- Fits phlebotomy trays
- Dual openings; vertical drop
- Dimensions: 6 <sup>7</sup>/<sub>16</sub>" H x 4 <sup>1</sup>/<sub>4</sub>" L x 4 <sup>1</sup>/<sub>4</sub>" W
- *Uses holder #400 to secure container on tray or tabletop*



**400 000**

**PHLEBOTOMY HOLDER**

**20/CASE**

- Translucent White allows visibility of container label
- Adhesive on bottom secures container to tabletop or tray
- Dimensions: 6 <sup>5</sup>/<sub>16</sub>" H x 4 <sup>1</sup>/<sub>4</sub>" L x 4 <sup>1</sup>/<sub>4</sub>" W
- *For use with phlebotomy container #100*



 **102 020**  
 **102 030**  
 **102 040**

**2-GALLON SHARPS CONTAINER**

**30/CASE**




**Translucent Beige**

**Translucent Red**

**Translucent Yellow**

- Accommodates syringes up to 60cc
- Built-in handles for easy transport
- Dimensions: 8 <sup>5</sup>/<sub>8</sub>" H x 11 <sup>5</sup>/<sub>8</sub>" L x 7 <sup>3</sup>/<sub>4</sub>" W
- *Uses bracket/key #435 and 440 for wall mounting*



 **303 020**  
 **303 030**  
 **303 040**

**3-GALLON SHARPS CONTAINER**

**12/CASE**

**Translucent Beige**

**Translucent Red**

**Translucent Yellow**

- Large capacity accommodates syringes up to 60cc
- Built-in handles for easy transport
- Dimensions: 13 <sup>1</sup>/<sub>2</sub>" H x 13 <sup>7</sup>/<sub>8</sub>" L x 6 <sup>7</sup>/<sub>8</sub>" W
- *Uses bracket #423 or 425 for wall mounting*