

VAC-U-STATION™



BEMIS
Fluid Management System

VAC-U-STATION™

For use with these VAC-U-PORT™ canisters

120ZD5P050	1200 cc Suction Canister
200ZD5P050	2000 cc Suction Canister
300ZD5P050	3000 cc Suction Canister

BEMIS
Fluid Management System

Bemis Manufacturing Corporation • 2500 East Main Street • Sheboygan Falls, WI 53085
USA: 920-452-4000 • 800-554-7800 • Fax: 920-452-6074 • Canada: 905-726-0000 x 4400

**A safe and economical way to dispose of suction canister contents.
An infectious waste management system for all areas of a healthcare
institution.**

VAC-U-STATION™ Fluid Management System

Dispose of suction canister contents directly into sanitary sewer system.




Vac-U-Port™ Suction Canisters

1200cc for Patient Rooms and Anesthesia
 2000cc for Operating Rooms and other areas
 3000cc for Operating Rooms

Canisters contain an Aerostat™ bacterial filter and float-valve shutoff. A complete line of accessories is also available.

Bemis VAC-U-STATION is used only in conjunction with Bemis Vac-U-Port Suction Canisters.



- * **Rapid Drain (250cc per second or more)**
- * **"ARM" to secure canister when draining**
- * **Secure disposal into a sanitary sewer system. No exposure (i.e. spilling, splashing, aerosolization) to healthcare professionals**
- * **Spray hose to rinse canister**
- * **User friendly handle is easy to activate**
- * **Viewport provides assurance that disposal channel is clear and clean**

- Eliminates the high cost of liquid infectious waste disposal, including expensive solidifiers.
- Small footprint - only 47.1"H x 18.7"W x 10.8"D needed to install.
- No electricity required. Attaches easily to existing plumbing.
- Usually mounted in soiled utility room.
- Daily maintenance requirements are minimal - approximately 3 minutes.

DESCRIPTION	CATALOG NO.	PACKAGING
Vac-U-Station™ Unit	VAC2000	1/case
1200cc Vac-U-Port Canister	1202DSP 050	48/case
2000cc Vac-U-Port Canister	2002DSP 050	48/case
3000cc Vac-U-Port Canister	3002DSP 050	36/case
Seal Maintenance Kit	VACSEAL KIT	1/each



Bemis Manufacturing Company • Health Care Products • 300 Mill Street, Sheboygan Falls, WI 53085
 920-467-4621 • 800-558-7651 • FAX 920-467-8573 • e-mail HCG@BemisMfg.com • www.BemisMfg.com