



Model 533910/534910
4-Canister Hi-Lo Stand

Case Pack: 1

Components:

Stand: Chrome finish steel – HRPO 1/16”
Coating Type: Nickel chrome finish
Plating thickness: 5ml

Wall Plate (4 each): Bayer 433 Lustran ABS

Manifold clip (2 each): Nylon

2” Casters (5 each): Steel and rubber

Height adjustment knob (1 each): Polypropylene

Maximum Height: 36” from top of stand to floor

Minimum Height: 26” from top of stand to floor



533910



534910

Note: Four of the 532710 plastic ring brackets are included with the 533910. These are not included with the 534910 as this model is designed to be used with the Quick-Fit liners and does not require separate brackets.