

healthVIE.comTM

leading the fight in infection prevention

Volume 1 / Issue 3

March 2010

A LOOK INSIDE

- hospital cleaning
- education modules
- disinfection practices
- prevention and control
- alcohol-based hand sanitizers
- liquid infectious waste disposal
- central sterile management

Reaching The Peak

of perioperative practice

healthVIE.com
P.O. Box 25310
Scottsdale, AZ 85255-0167

Sharps Safety

- BD Multi-use nestables; BD, www.bd.com.
- DIALYSAFETY™ Dialysis Sharps Disposal Containers; Covidien, www.kendallsharpsafety.com.

Sharps Containers, Disposable

Numerous areas in healthcare need to have a sharps collectors that meet their specific needs. The most commonly used sharps containers are disposable in nature; they should be disposed of when they are three-quarters full.

- BD Next Generation In-Room Sharps Collector; BD, www.bd.com.
- Complete In-room System; Health Care Logistics Inc., www.healthcarelogistics.com.
- EPS Sharps Disposal; Medidose Inc., www.medidose.com.
- GatorGuard Sharps Disposal Containers, SharpStar Sharps Disposal Containers, RENEWABLES™ Sharps Disposal Containers manufactured using recycled plastics, GATORGUARD™ In-Patient Room Sharps Containers, IN-ROOM™ System with Always-Open Lids, IN-ROOM System with Mailbox-Style Lids, IN-ROOM System with SHARPSTAR™ Counterbalance Door, IN-ROOM System with SHARPSTAR Lift Assisted Lid, SHARPS-A-GATOR™ Tortuous Path Sharp Containers, SHARPSTAR IV 3 & 4 Gallon Sharps Container with SHARPSTAR Counterbalance Door, Horizontal-Drop Opening Lids, MONOJECT™ Sharps Containers, Multi-Purpose Sharps Container, SHARPS-A-GATOR Sharps Collection & Disposal Systems, ENVIROSTAR™ non-infectious waste sharps disposal container; Covidien, www.kendallsharpsafety.com.
- Multi-Safe Disposal Box; Sarstedt, www.sarstedt.com.
- Sharps Disposal Containers; Greiner Bio-One, www.gbo.com/preanalytics.
- Sharps Secure™ sharps containers, Pro-Tec® sharps containers; Sharps Compliance Inc., www.sharpsinc.com.
- Sharps Sentinel Cover, Patient Room System, Phlebotomy & Multi-Use Containers, Wall Cabinet System, Large Volume and Chemotherapy Containers; Bemis Medical, www.bemis.com.

Sharps Containers, Large Volume

Healthcare workers and others need to be able to use a sharps disposal container to place needles and other sharps into as mandated by law. Large products such as endoscopic instruments need a safe place for disposal.

- Large Sharps Containers; BD, www.bd.com.
- Large Volume Sharps Containers; Bemis, www.Bemis.com.
- Large Volume Sharps Container; Health Care Logistics, www.healthcarelogistics.com.
- Large Volume Containers with Hinged, Sealing Gasket Lids, Large Volume Containers with Sliding, Sealing

Gasket Lids, SHARPS-A-GATOR™ Large Volume Sharps Containers; Covidien, www.kendallsharpsafety.com.

- Sharps Disposal Containers; Greiner Bio-One, www.gbo.com/preanalytics.



Bemis large volume sharps containers, designed for high volume areas, can easily accommodate large sharps, needle counters, trocars, endoscopic instruments, etc. Two different lid styles allow for temporary closure when not in use and one-handed disposal by activating the foot pedal.

Sharps Containers, Mail-back

Mail-back sharps containers help to provide a safe, cost-effective method of containment and disposal of medical waste for many light users of sharps including home healthcare, dental and physician offices, satellite healthcare facilities, health clubs, and hotels. Simply put the sharps into the container, seal it up, and ship it! More than eight million people