

# CONTOUR™ RIGID LIPOSUCTION CANISTERS

DESIGNED FOR CONFIDENT USE DURING LIPOSUCTION PROCEDURES



## FEATURES + BENEFITS

Contour is designed for confident use during liposuction and other procedures requiring high vacuum pressure and/or the presence of body lipids. For an airtight seal during procedure, Contour features a self-sealing lid with port connection flexibility that can be used with or without elbow connections.

- Manufactured with lipid-resistant materials to prevent the microscopic breakdown of plastic that cause canister malfunctions
- Canister offers excellent clarity for enhanced visualization of contents with easy to read graduations
- Hi-Flow lid design provides high flow rate with mechanical shut off valve for prevention of overflow and vacuum line contamination
- The integrated Aerostat filter provides >99.99% bacterial efficiency\*

Size	Lipid Resistant	Peak Vacuum Pressure
2000cc	•	29inHg / 736mmHg
3000cc	•	29inHg / 736mmHg

Disposable, single use, non-sterile  
Not made with natural rubber latex; BPA and DEHP free  
Includes two elbows with Contour canisters

## CONTOUR SHUTOFF + FILTER



- A Aerostat Filter** provides >99.99% bacterial efficiency\*
- B Filter Retainer**
- C Mechanical Shutoff Valve** with soft gasket seal
- D Rigid Float Support** keeps float valve securely in place

## 2000 + 3000CC CONTOUR LID



- E Tandem Port** for multiple canister connections
- F Patient Port** stepped design accommodates various tubing sizes or can be used with elbow to prevent kinking
- G Port Caps** attached to all lids for quick and secure capping
- H Vacuum Port** for use with or without elbow
- I Accessory Port**

\* Data on file

## CONTOUR LIPOSUCTION CANISTERS - PINK LIDS

PRODUCT/DESCRIPTION	CATALOG#	PACKAGING
 <p><b>2000CC CONTOUR LIPOSUCTION CANISTER</b></p> <ul style="list-style-type: none"><li>• Float valve shutoff with Aerostat filter, two elbows</li></ul> <p><b>Bottom Diameter: 4.125in Top Diameter: 6.25in Height: 9.875in</b></p>	2000C 025	12/CASE
 <p><b>3000CC CONTOUR LIPOSUCTION CANISTER</b></p> <ul style="list-style-type: none"><li>• Float valve shutoff with Aerostat filter, two elbows</li></ul> <p><b>Bottom Diameter: 4.375in Top Diameter: 6.25in Height: 12.875in</b></p>	3000C 025	12/CASE

### Bemis Manufacturing Company

Health Care Division | 300 Mill Street | PO Box 901 | Sheboygan Falls, WI 53085 USA  
Email: [Customer.Service@BemisHealthCare.com](mailto:Customer.Service@BemisHealthCare.com) | Phone: 920.467.4621 | Fax: 800.499.8160

4/2025 Rev 02  
©2025 0B7752121

