

# CONTAINER MOUNTING OPTIONS QUICK REFERENCE TABLE

Sentinel™  
by BEMIS®

CONTAINER CATALOG #	COLOR	SIZE	MOUNTING TYPE & CATALOG #						
			5-QT SHARPS CABINET	3-GAL SHARPS CABINET	PHLEBOTOMY HOLDER	WIRE BRACKET (NON-LOCKING)	WIRE BRACKET (LOCKING)*	WALL BRACKET & KEY**	ROLLING CART
			1051	1031	400 000	423 000	425 000	445 020	430 030
<b>SHARPS CONTAINERS</b>									
100 040	Yellow	1-Quart			•				
100 030	Red	1-Quart			•				
175 040	Yellow	5-Quart	•					•	
175 030	Red	5-Quart	•					•	
175 020	Beige	5-Quart	•					•	
102 040	Yellow	2-Gallon						•	
102 030	Red	2-Gallon						•	
102 020	Beige	2-Gallon						•	
303 040	Yellow	3-Gallon				•	•		
303 030	Red	3-Gallon				•	•		
303 020	Beige	3-Gallon				•	•		
333 040	Yellow	3-Gallon		•		•	•		
333 030	Red	3-Gallon		•		•	•		
333 020	Beige	3-Gallon		•		•	•		
108 030	Red	8-Gallon							•
108 040	Yellow	8-Gallon							•
111 030	Red	11-Gallon							•
111 040	Yellow	11-Gallon							•
<b>CHEMOTHERAPY CONTAINERS</b>									
202 004	White	2-Gallon						•	
202 040	Yellow	2-Gallon						•	
208 004	White	8-Gallon							•
208 040	Yellow	8-Gallon							•
211 004	White	11-Gallon							•
211 040	Yellow	11-Gallon							•
<b>NON-HAZARDOUS PHARMACEUTICAL CONTAINERS</b>									
4002 050	Blue	2-Gallon						•	
4008 050	Blue	8-Gallon							•
4011 050	Blue	11-Gallon							•
<b>HAZARDOUS PHARMACEUTICAL CONTAINERS (RCRA)</b>									
5002 070	Black	2-Gallon						•	
5008 070	Black	8-Gallon							•
5011 070	Black	11-Gallon							•

\* Purchase padlock and key separately - catalog # 426 000

\*\* Wall brackets and keys can be purchased separately - catalog # 435 020 and 440 020

## Bemis Manufacturing Company

Health Care Division | 300 Mill Street | P.O. Box 901 | Sheboygan Falls, Wisconsin 53085 USA

Email: Customer.Service@BemisHealthCare.com | Phone: 920.467.4621 | Toll-Free Phone: 800.558.7651 | Fax: 800.499.8160

3/2022 Rev 01

©2022 0B7756011

