

BEMIS
HEALTH CARE™



CONTOUR^{BY BEMIS} LipoSuction Canister™

THE DIFFERENCE IS
CLEAR

CONTOUR™ LIPOSUCTION CANISTERS

FEATURES

- Engineered for high vacuum liposuction procedures
- Designed and manufactured using lipid-resistant material
- Excellent clarity for enhanced visualization and accuracy
- Large, easy to read graduations
- Lids contain multiple ports for versatility and tandem set up
- High performance filter and shutoff protects against contamination and overflow
- Available in 2 sizes to accommodate varying needs
- Nestable for storage convenience
- Disposable, single use
- May be used in conjunction with Bemis Quick-Drain™
- Vacuum level rating of 29" Hg/736mm Hg

PRODUCT/DESCRIPTION

CATALOG

PACKAGING

CONTOUR™ LIPOSUCTION CANISTERS



2000CC CONTOUR™ LIPOSUCTION CANISTER

with Aerostat® filter, float valve shutoff

Bottom Diameter: 4" **Top Diameter:** 6 1/8" **Height:** 9 3/4"

2000C 025

12/CASE



3000CC CONTOUR™ LIPOSUCTION CANISTER

with Aerostat® filter, float valve shutoff

Bottom Diameter: 4" **Top Diameter:** 6 1/8" **Height:** 12 3/4"

3000C 025

12/CASE