

QUICK-FIT™ SUCTION LINERS

For facilities that prefer the convenience of a liner system



FEATURES

Aerostat® filter, Hi-flow shutoff system with anti-reflux valve

Latex-free


Available in two sizes, 1500cc and 3000cc

Liner comes individually bagged and compressed to minimize storage space

Pre-assembled sealed lids and bags save time and eliminate confusion

Reusable outer canister is made of durable polycarbonate

EASY TO USE



- A** Color-coded lids and brackets allow easy visual selection.
- B** Innovative rib design allows easy insertion and removal.
- C** Elbow allows quick hookup.
- D** Available with or without built-in stopcock.

- E** Stepped patient port permits direct patient tubing connection.
- F** Large port for attachment of Specimen Collector when necessary.
- G** Built-in handle facilitates removal.

PRODUCT/DESCRIPTION

CATALOG #

PACKAGING

QUICK-FIT™ SUCTION LINER SYSTEM

INDIVIDUALLY BAGGED SUCTION LINERS



1500CC SUCTION LINER

with attached yellow lid, Aerostat® filter, float valve shutoff, anti-reflux valve, one elbow

Bottom Diameter: 4 1/8" **Top Diameter:** 6 1/8" **Height:** 8 3/4"

1504

50/CASE



3000CC SUCTION LINER

with attached pink lid, Aerostat® filter, float valve shutoff, anti-reflux valve, one elbow

Bottom Diameter: 3 3/4" **Top Diameter:** 6 1/8" **Height:** 15 3/4"

3004

50/CASE

REUSABLE OUTER CANISTERS *Available with or without stopcocks*



1500CC REUSABLE OUTER CANISTER

with attached yellow bracket

Bottom Diameter: 4 5/8" **Top Diameter:** 6 1/8" **Height:** 7 1/2"

1500 10

10/CASE

1500 01

1/CASE



1500CC REUSABLE OUTER CANISTER WITH STOPCOCK

with attached yellow stopcock bracket

Bottom Diameter: 4 5/8" **Top Diameter:** 6 1/8" **Height:** 7 1/2"

1500SC 10

10/CASE

1500SC 01

1/CASE



3000CC REUSABLE OUTER CANISTER

with attached pink bracket

Bottom Diameter: 4 1/8" **Top Diameter:** 6 1/8" **Height:** 14 3/8"

3000 10

10/CASE

3000 01

1/CASE



3000CC REUSABLE OUTER CANISTER WITH STOPCOCK

with attached pink stopcock bracket

Bottom Diameter: 4 1/8" **Top Diameter:** 6 1/8" **Height:** 14 3/8"

3000SC 10

10/CASE

3000SC 01

1/CASE