

# QUICK-FIT™ SUCTION LINERS



For facilities that prefer the convenience of a liner system



## FEATURES & BENEFITS

Three bag liner sizes available and designed for strength with multilayer nylon and polyethylene materials

Transparent liner for clear visualization of contents

Integrated Aerostat® filter has a >99.99% bacterial efficiency of aerosolized particulate\*\*

Hi-Flow shut off valve closes automatically to prevent overflow and vacuum line contamination

Anti-reflux valve on patient port prevents fluid escape upon disconnection

Reusable outer canister made of durable polycarbonate with the canister bracket and liner lids color coded to quickly identify size

May be used in conjunction with Bemis Quick-Drain™

Compact liner design for space savings in storage

Latex, BPA, and DEHP<sup>1</sup> free

Liner: single use, disposable; outer canister: reusable

## EASY TO USE



**A** Color-coded lids and brackets allow easy visual selection.

**B** Strong nylon/polyethylene waste bag.

**C** Patient Port elbow for prevention of kinked tubing.

**D** Stopcock for easy on/off of vacuum source.

**E** Patient port with Anti-Reflux valve - use with or without elbow.

**F** Large pour spout for easy emptying and for attachment of Specimen Adaptor when necessary.

**G** Built-in handle facilitates removal and transport.

**H** Tandem hookup port for multiple liner configuration.

**QUICK-FIT™ SUCTION LINER SYSTEM**

**BULK PACK SUCTION LINERS**



**1000CC SUCTION LINER, WHITE LID**  
 Aerostat® filter, float valve shutoff, patient port  
 anti-reflux valve, one elbow  
**Bottom Width: 7" Top Diameter: 6 3/8" Height: 7 7/8"**

1004RBM N 50/CASE



**1500CC SUCTION LINER, YELLOW LID**  
 Aerostat® filter, float valve shutoff, patient port  
 anti-reflux valve, one elbow  
**Bottom Width: 7" Top Diameter: 6 3/8" Height: 10 1/4"**

1504RBM N 50/CASE



**3000CC SUCTION LINER, PINK LID**  
 Aerostat® filter, float valve shutoff, patient port  
 anti-reflux valve, one elbow  
**Bottom Width: 6 7/8" Top Diameter: 6 3/8" Height: 17 5/8"**

3004RBM N 50/CASE

**REUSABLE OUTER CANISTERS**



**1000CC REUSABLE OUTER CANISTER<sup>1</sup>**  
 with attached white stopcock bracket  
**Bottom Diameter: 4 3/4" Top Diameter: 6" Height: 5"**

1000SC2 10 10/CASE  
 1000SC2 01 1/CASE



**1500CC REUSABLE OUTER CANISTER WITH STOPCOCK**  
 with attached yellow stopcock bracket  
**Bottom Diameter: 4 5/8" Top Diameter: 6" Height: 7 3/8"**

1500SC 10 10/CASE  
 1500SC 01 1/CASE



**3000CC REUSABLE OUTER CANISTER WITH STOPCOCK**  
 with attached pink stopcock bracket  
**Bottom Diameter: 4 1/2" Top Diameter: 6" Height: 14 1/4"**

3000SC 10 10/CASE  
 3000SC 01 1/CASE

<sup>1</sup> 1000CC Reusable Outer Canister connector hose contains presence of phthalate (DEHP)