

BEMIS SUCTION CANISTER AND LINER WEIGHTS & DIMENSIONS

SIZE/DESCRIPTION	BOTTOM DIAMETER (IN)	TOP DIAMETER (IN)	HEIGHT (IN)	EMPTY WEIGHT (LBS)	FULL WEIGHT (LBS)
US CUSTOMARY					
800CC Hydrophobic	3 ½	5 ⅜	7 ¼	0.32	2.08
800CC Hi-Flow	3 ½	5 ⅜	6 ⅝	0.33	2.09
1200CC Rigid**	4 ⅞	6 ¼	7 ½	0.53	3.17
2000CC Rigid**	4 ⅞	6 ¼	9 ⅞	0.66	5.07
3000CC Rigid**	4 ⅞	6 ¼	12 ⅞	0.86	7.47
1000CC Liner	7 (width)	6 ⅜	7 ⅞	0.27	2.47
1500CC Liner	7 (width)	6 ⅜	10 ¼	0.28	3.59
3000CC Liner	6 ⅞ (width)	6 ⅜	17 ⅝	0.29	6.90
1000CC Reusable Outer Canister*	4 ¾	6	5	N/A	N/A
1500CC Reusable Outer Canister*	4 ⅞	6	7 ⅞	N/A	N/A
3000CC Reusable Outer Canister*	4 ⅞	6	14 ¼	N/A	N/A
*NOTE: The bracket attached to the Quick-Fit Reusable Outer Canister is 2 ⅝" long					
**NOTE: "Rigid" signifies all rigid canisters that contain listed volumes e.g., Hi-Flow, Hydrophobic, Contour					

SIZE/DESCRIPTION	BOTTOM DIAMETER (MM)	TOP DIAMETER (MM)	HEIGHT (MM)	EMPTY WEIGHT (KG)	FULL WEIGHT (KG)
METRIC					
800CC Hydrophobic	88.9	136.5	184.2	0.14	0.94
800CC Hi-Flow	88.9	136.5	169.5	0.15	0.95
1200CC Rigid**	104.8	155.6	190.5	0.24	1.44
2000CC Rigid**	104.8	155.6	250.8	0.30	2.30
3000CC Rigid**	111.1	155.6	327	0.39	3.39
1000CC Liner	177.8 (width)	161.9	200	0.12	1.12
1500CC Liner	177.8 (width)	161.9	260.4	0.13	1.63
3000CC Liner	174.6 (width)	161.9	447.7	0.13	3.13
1000CC Reusable Outer Canister*	120.7	152.4	127	N/A	N/A
1500CC Reusable Outer Canister*	117.5	152.4	187.3	N/A	N/A
3000CC Reusable Outer Canister*	104.8	152.4	362	N/A	N/A
*NOTE: The bracket attached to the Quick-Fit Reusable Outer Canister is 59.375 mm long					
**NOTE: "Rigid" signifies all rigid canisters that contain listed volumes e.g., Hi-Flow, Hydrophobic, Contour					