

## Bemis Suction Canister and Liner Weights and Dimensions

<u>Size/Description</u>	<u>Bottom Diameter (in.)</u>	<u>Top Diameter (in.)</u>	<u>Height (inches)</u>	<u>Empty Wt. (lbs)</u>	<u>Full Wt. (lbs)</u>
800cc Hydrophobic	3-3/8"	5"	7-1/4"	0.39	2.09
800cc Hi-Flow	3-3/8"	5"	6-5/8"	0.39	2.09
1200cc Rigid	4"	6-1/8"	7-1/2"	0.77	3.37
1200cc Vac-U-Port	4"	6-1/8"	8"	0.77	3.37
2000cc Rigid	4"	6-1/8"	9-3/4"	0.79	5.19
2000cc Vac-U-Port	4"	6-1/8"	10-1/4"	0.79	5.19
3000cc Rigid	4-3/8"	6-1/8"	12-3/4"	0.87	7.47
3000cc Vac-U-Port	4-3/8"	6-1/8"	13-1/2"	0.87	7.47
1500cc Liner	4-1/4"	6-1/8"	8-3/4"	0.28	3.9
3000cc Liner	3-3/4"	6-1/8"	15-3/4"	0.32	7.1
1500cc Reusable Outer Canister*	4-5/8"	6-1/8"	7-1/2"	N/A	N/A
3000cc Reusable Outer Canister*	4-1/8"	6-1/8"	14-3/8"	N/A	N/A

\*NOTE: The bracket attached to the Reusable Outer Canister for Liners is 2-3/8" long

### METRIC MEASUREMENTS

<u>Size/Description</u>	<u>Bottom Diameter (mm)</u>	<u>Top Diameter (mm)</u>	<u>Height (mm)</u>	<u>Empty Wt. (kg)</u>	<u>Full Wt. (kg)</u>
800ml Hydrophobic	84.375	125	181.25	0.1755	0.9405
800ml Hi-Flow	84.375	125	165.625	0.1755	0.9405
1200ml Rigid	100	153.125	187.5	0.3465	1.5165
1200ml Vac-U-Port	100	153.125	200	0.3465	1.5165
2000ml Rigid	100	153.125	243.75	0.3555	2.3355
2000ml Vac-U-Port	100	153.125	256.25	0.3555	2.3355
3000ml Rigid	109.375	153.125	318.75	0.3915	3.3615
3000ml Vac-U-Port	109.375	153.125	337.5	0.3915	3.3615
1500ml Liner	106.25	153.125	218.75	0.126	1.755
3000ml Liner	93.75	153.125	393.75	0.126	3.195
1500ml Reusable Outer Canister*	115.625	153.125	187.5	N/A	N/A
3000ml Reusable Outer Canister*	103.125	153.125	359.375	N/A	N/A

\*NOTE: The bracket attached to the Reusable Outer Canister for Liners is 59.375mm long