



FEATURES

Gasketed lids and absorbent spill pads provide leak resistance.

Available in three sizes.

Clearly labeled with chemotherapy warnings.

Nestable to save space and reduce shipping costs.

Multilingual labeling ensures user safety.

Note: Chemotherapy containers must be incinerated.

CHEMOTHERAPY SHARPS CONTAINERS

Designed for safe collection and containment of trace chemotherapy waste including needles, syringes, IVs, empty or partially empty vials, as well as gloves and gowns. The EPA defines "trace" amounts of waste as containing no more than 3% by weight of the original contents.





■ 202 004
■ 202 040

2-GALLON CHEMOTHERAPY CONTAINER 30/CASE

White w/gasket lid/absorbent pad
Yellow w/gasket lid/absorbent pad

- Screw-on cap; Posts on lid for holding cap in place when in use
- For use on tabletop or wall-mounted
- Dimensions: 9" H x 11 5/8" L x 7 3/4" W
- *Uses bracket/key #435 and 440 for wall mounting*



■ 208 004
■ 208 040

8-GALLON CHEMOTHERAPY CONTAINER 10/CASE

White w/gasket lid/absorbent pads
Yellow w/gasket lid/absorbent pads

- Separate round opening accepts small items
- Hinged lid accommodates large items
- Molded-in handles for easy lifting and transport
- Dimensions: 15 7/8" H x 16 1/2" L x 11 13/16" W
- *Uses cart #430 for mobility*



■ 211 004
■ 211 040

11-GALLON CHEMOTHERAPY CONTAINER 6/CASE

White w/gasket lid/absorbent pads
Yellow w/gasket lid/absorbent pads

- Separate round opening accepts small items
- Hinged lid accommodates large items
- Molded-in handles for easy lifting and transport
- Dimensions: 22 1/2" H x 16 1/2" L x 11 13/16" W
- *Uses cart #430 for mobility*